

MEDICAL PACKAGING

BOLSAPACK • PACKAGING FOR LOW TEMPERATURE PLASMA (H₂O₂)/EO

Pouches and rolls, manufactured with Tyvek® and plastic film (transparent). Indicator for low temperature is printed in packaging.

TYVEK®/FILM FLAT BAG


Code	Sizes		U/Box	Box/cm	Box/Kg
MBT1160	75 x 210 mm	3,0" x 8,3"	1.000	16 x 22 x 15	2,4
MBT1161	100 x 280 mm	3,9" x 11,0"	1.000	16 x 22 x 15	4,1
MBT1162	150 x 320 mm	5,9" x 12,6"	1.000	31 x 33 x 15	7
MBT1163	200 x 420 mm	7,9" x 16,5"	350	43 x 21 x 11	4
MBT1164	250 x 500 mm	9,8" x 19,7"	350	28 x 60 x 12	6,5
MBT1165	320 x 560 mm	12,6" x 22,0"	400	41 x 59 x 11	10,5
MBT1166	400 x 580 mm	15,7" x 22,8"	350	41 x 59 x 11	12

TYVEK®/FILM FLAT ROLL

Code	Sizes		U/Box	Box/cm	Box/Kg
MRT1171	75 mm x 70 m	3,0" x 230 Feet	4	31 x 25 x 22,5	3,5
MRT1172	100 mm x 70 m	3,9" x 230 Feet	4	41 x 21 x 22	4,5
MRT1174	150 mm x 70 m	5,9" x 230 Feet	4	37 x 34,0 x 21,5	6,5
MRT1175	200 mm x 70 m	7,9" x 230 Feet	2	41 x 21 x 22	4,5
MRT1176	250 mm x 70 m	9,8" x 230 Feet	2	43 x 26 x 22,5	6,0
MRT1177	300 mm x 70 m	11,8" x 230 Feet	2	37 x 21 x 34	6,0
MRT1180	350 mm x 70 m	13,8" x 230 Feet	1	41 x 21,8 x 23	3,6
MRT1178	400 mm x 70 m	15,7" x 230 Feet	1	41 x 21 x 22	4,5
MRT1179	500 mm x 70 m	19,7" x 230 Feet	1	52,5 x 18 x 18	6,0



Tyvek®/Film bag



Opening instructions



Tyvek®/Film roll



Plasma indicator Tyvek®

INDUSTRIAL PACKAGING

Although our production line is focused in supplying stock products, we also print medical paper and films for manufacturers of medical devices for their packaging machines. Our experience and knowledge in manufacturing stock items, is an added value for our industrial customers, providing exactly what customer needs.

PRODUCTS

- Medical paper
- Pa/pp, pet/pp, pet/pe, pa/pe films
- Tyvek®
- Private label

(Tyvek® is a registered trade mark of DUPONT)